## Infrarotleuchte IRL

### Infrared lamp IRL

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Lampada a raggi infrarossi IRL

Lâmpara infrarroja IRL

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Infraroodlamp IRL

Infrapunalamppu IRL

Infraröd lampa IRL

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### 1.1 Thank you very much

Thank you for your confidence in us and congratulations!

The **IRL** infrared lamp that you have purchased is a top quality product from **MEDISANA**.

This device is for providing cosmetic and therapeutic treatment using extremely intensive infrared heat that penetrates deep into the skin.

In order to achieve the desired effect with your **IRL** infrared lamp from **MEDISANA** in the long term, we recommend that you read the following information on its use and maintenance carefully.

**Read the instruction manual carefully before using this device, especially the safety instructions, and keep the instruction manual for future use.**

**Should you give this device to another person, it is vital that you also pass on these instructions for use.**

### 1.2 Important

- Use the device only according to its intended purpose as specified in the instruction manual.
- The warranty will be invalidated if the device is used for purposes other than those for which it is intended.
- Do not use the device if it is not working correctly, if it has been dropped or has fallen into water or has been damaged.
- Do not use the device if the cable or socket is damaged and keep the cable away from hot surfaces.
- Do not attempt to repair the device yourself in the event of a malfunction since this will invalidate the warranty. Repairs should only be carried out by authorised service centres.
- The device is not determined for commercial use.
- Before using make sure that the connecting cable is not damaged. It must not come into contact with the hot electric light bulb.
- Before connecting the device to your power supply, please ensure that the supply voltage stated on the rating plate is compatible with your mains supply.
- Place device on a firm, level surface during use so that it cannot tip over. Do not attach device to walls or ceilings.
- Do not leave the device unattended when it is switched on.
- Keep flammable objects such as textiles, glasses and combs away from the infrared lamp.
- Keep the infrared lamp at least 80 cm away from flammable objects.
- Do not use the device unless it is completely dry.
- Do not use the device outdoors or in damp rooms.
- The device must not be used in the immediate vicinity of bathtubs, showers or washbasins filled with water, and must not be used with wet hands.
- Protect the device from moisture.
- Do not immerse the device in water or other fluids.
- If, nonetheless, liquid should penetrate the device, unplug the mains cable immediately.

## 1.2 Important

- Always ensure that your hands are dry when you get hold of the lamp and/or the plug.
- Never carry, pull or rotate the unit by the mains cable and never pinch the cable between other objects.
- Do not cover the device when it is warm.
- Protect the device from shocks and vibration.
- **Attention, scalding danger!**
  - The bulb becomes hot.
  - Do not touch the infrared bulb and do not remove with the power cord plugged in.
  - Always disconnect the device from the power supply and allow the lamp to cool before touching it.
- Always pull on the plug itself when disconnecting the device from the power supply. Avoid pulling on the power cord!
- Always switch the device off after use, before cleaning and care and in the event of malfunctions during operation by disconnecting the plug from the socket.
- Allow the infrared lamp to cool before putting it away.
- Children may not use the device. Keep the device away from children. Medical products are not designed to be used as a toy!
- If the unit is to be used by or on people who are handicapped in any way or in the presence of children, this must always be done under supervision.
- If the scale is handed to third parties, please be sure to give them this instruction manual.

## 1.3 Important advice - for your health's sake

- Treatment with the **MEDISANA** infrared lamp **IRL** cannot replace diagnoses and treatment by a physician.
- Treatment after injuries or operations may only be given under the supervision and guidance of a doctor.
- Pain-killers and relaxants reduce the sensitivity of the skin to pain and heat; the same applies to the consumption of alcohol. Consult your doctor before using the device.
- Persons, whose sensitivity is disturbed or impaired may only apply the device according to the instructions of a physician.
- This device is unsuitable for persons who are insensitive to heat.
- Do not treat any parts of the body that have swellings, burns, inflammation, peeling skin, wounds or other sensitive points.
- Do not use the device on skin that feels painful or numb.
- If the skin becomes irritated, discontinue the treatment and consult a physician.
- Always seek the advice of a doctor in case of doubt and follow his instructions for treatment.
- The distance to the lamp must be at least 40 cm to avoid that the intensity of radiation gets too high and the warmth is felt as painful.

- Always close your eyes during radiation in order to prevent retinal injuries; do not look directly into the infrared light. You can also cover your eyes with a swab.
- Avoid using the lamp for long periods, since this can result in burning.
- Do not exceed the maximum radiation time of 15 minutes.
- The treatment should be pleasant. Stop at once and consult your doctor if you feel pain or if the treatment seems unpleasant.

#### 1.4 Scope of delivery and packaging

Please check first whether the device is complete.

The following belong to the scope of delivery:

- 1 **MEDISANA** infrared lamp **IRL**
- 1 instruction manual

The packaging can be reused or recycled.

Please dispose properly of any packaging material no longer required.

If you notice any transport damage during unpacking, please contact your dealer without delay.

#### **WARNING**

**To avoid risk of suffocation, keep packaging film away from children.**

#### 2.1 How does infra-red radiation work?

Heat promotes blood circulation and supports the healing process. The infra-red lamp generates intensive radiation that is led to the part of the body that is being treated in the form of light. This light beam has the property of penetrating down to the deepest skin layers. The radiation is then converted into intensive, extremely effective heat in the tissue and muscles. The heat that is generated in this way expands the vessels, encourages blood circulation and stimulates the circulation and metabolism. The pores in the skin also expand. Infra-red radiation is therefore suitable for both therapeutic and cosmetic treatment. The stimulation of the metabolism reinforces the body's internal defences and accelerates healing processes.

### 3.1 Application area of the MEDISANA IRL infrared lamp

#### ***Therapeutic treatment***

The **MEDISANA IRL** infrared lamp supports domestic therapy and can accompany medical treatment given by a doctor, such as:

- Colds (consult a doctor if you have a high temperature).
- Tension
- Cramp
- Muscular pain
- Rheumatic complaints
- Treatment of ear, nose and throat complaints

#### ***Cosmetic treatment***

- The **MEDISANA IRL** infrared lamp is suitable for assisting facial and beauty care, particularly in the case of bad skin.
  - Remove all traces of make-up before treatment.
  - The effect of the infrared radiation opens up the pores in the skin.
  - This results in better and more effective cleaning of the skin.
  - Creams and other cosmetic materials are absorbed more intensively and quickly by the skin.

No side-effects are known for the treatments mentioned in these usage instructions, provided that all of the safety instructions are followed (see pages 8 – 10). If you do not see an improvement after several treatments, please contact a doctor.

In all other cases that are not mentioned here, infrared radiation therapy must not be carried out before consulting a doctor.

### 3.2 Operation

- Please follow the safety instructions.
- Check the connecting cable **2**.
- Place the device on its base **3** on a firm, level surface.
- The infrared bulb **4** must be firmly screwed in.
- The device has an adjustable lamp housing **1**. The 4 settings from 0 ° to 45 ° allow the device to be individually positioned in accordance with requirements.
- In order to prevent burns, always grip the lamp housing **1** from behind if you wish to change the position.
- Plug the device into the mains and it will start up immediately. Please operate the equipment at a distance at which the heat feels comfortable, but no less than 40 cm. Increase the distance if the heat feels uncomfortable.
- The duration of use is oriented to your sensitivity, but do not use for more than 15 minutes.
- Unplug the equipment from the mains when treatment is complete.

**IMPORTANT NOTES**

The duration of use should be 5 to 10 minutes, once or twice per day. Do not exceed the maximum radiation time of 15 minutes per session under any circumstances.

The infrared radiation emitted by the device is not suitable for tanning.

#### 4.1 Cleaning and maintenance

- Before cleaning and replacing the bulb of the infrared lamp, unplug the equipment from the mains and allow the device to cool.
- Clean the device after cooling, ① and ③, with a damp cloth and then rub dry.
- Do not use aggressive cleaning and scouring materials.
- Never immerse the device in water.
- Please ensure that water does not penetrate the internals of the device. Do not use the device again until it is completely dry.
- Always replace the bulb of the infra-red lamp with the same bulb type.
- Store the device in a dry, cool place.

#### *infrared light bulb* ④

It must not be used one with a higher power as indicated on the rating plate. Replacements for the light bulb can be bought where the appliance has been bought or ordered directly from the producer (indicate model!).

The following replacements for the light bulb have to be used:  
GE PAR38 IR 150 W or Philips Infrared PAR38E 150 W

In case the glass of the light bulb is broken immediately disconnect the appliance. Attention! Injuring danger due to splinters of glass. Let the appliance cool down. Carefully remove the part of the light bulb which is screwed in the socket with a pair of pincers (protect your eyes) or have the appliance repaired by a service agent.

The light bulb is not part of the guarantee, it won't be replaced free of charge.

#### 4.2 Disposal

This product must not be disposed together with the domestic waste.

All users are obliged to hand in all electrical or electronic devices, regardless of whether or not they contain toxic substances, at a municipal or commercial collection point so that they can be disposed of in an environmentally acceptable manner.

Consult your municipal authority or your dealer for information about disposal.

**4.3****Technical Data**

Name und model	:	<b>MEDISANA</b> infrared lamp IRL
Current supply	:	230 V~ / 50 Hz
Consumption	:	150 Watt
Dimensions L x W x H	:	approx. 16 x 14 x 20 cm
Weight	:	approx. 0,6 kg
Safety class	:	II 
Article no.	:	88252
EAN-Number	:	40 15588 88252 4
Spare bulb	:	88255

**CE 0297**

**In accordance with our policy of continual product improvement, we reserve the right to technical and design changes without notice.**

**Electromagnetic compatibility:**

The device complies with the EN 60601-1-2 (2001) standard for electromagnetic compatibility.

Inquire at **MEDISANA** for details on this measurement data.

### 5.1 Warranty and repair terms

Please contact your dealer or the service centre in case of a claim under the warranty. If you have to send in the unit, please enclose a copy of your receipt and state what the defect is.

The following warranty terms apply:

1. The warranty period for **MEDISANA** products is three years from date of purchase. In case of a warranty claim, the date of purchase has to be proven by means of the sales receipt or invoice.
2. Defects in material or workmanship will be removed free of charge within the warranty period.
3. Repairs under warranty do not extend the warranty period either for the unit or for the replacement parts.
4. The following is excluded under the warranty:
  - a. All damage which has arisen due to improper treatment, e.g. non-observance of the user instructions.
  - b. All damage which is due to repairs or tampering by the customer or unauthorised third parties.
  - c. Damage which has arisen during transport from the manufacturer to the consumer or during transport to the service centre.
  - d. Accessories which are subject to normal wear and tear.
5. Liability for direct or indirect consequential losses caused by the unit are excluded even if the damage to the unit is accepted as a warranty claim.

#### **MEDISANA AG**

Itterpark 7-9  
40724 Hilden  
Germany

Tel.: +49 (0) 2103 / 2007-60  
Fax: +49 (0) 2103 / 2007-626  
eMail: [info@medisana.de](mailto:info@medisana.de)  
Internet: [www.medisana.de](http://www.medisana.de)

The service centre address is shown on the last page.

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D&L Marketing Oy
ALERE Electronics AB