Personenwaage PST Personal scales PST

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Bilancia pesapersona PST

Báscula de baño PST

Balança PST

Personenweegschaal PST

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Personvåg PST

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1.1 Thank you very much

Thank you for your confidence in us and congratulations!

By purchasing the personal scales **PST**, you have acquired a quality product from **MEDISANA**.

This unit is intended for weighing people.

To ensure that your **MEDISANA** personal scales **PST** provide successful results for a long time, we recommend that you read carefully through the following notes on their use and care.

1.2 Important

Please read the following safety information carefully before using the device and keep these instructions in a safe place for later use.

- Any treatment or diet for people who are underweight or overweight requires the qualified advice of a professional (doctor, dietician). The values determined with the scales may be source of useful support.

- ***Risk of tipping over!***

Do not place the scale on an uneven surface.

Never step onto a corner of the scales.

When weighing yourself, do not stand on one side or the edge of the scales.

- ***Risk of slipping!***

Do not stand on the scale with wet feet.

Do not stand on the scale wearing socks.

- Use the device only according to its intended purpose as specified in the instruction manual.
- The warranty will be invalidated if the device is used for purposes other than those for which it is intended.
- The scale has been manufactured for domestic use. It is not suitable for commercial use in hospitals or other medical institutions.
- Do not use the device if it is not working correctly, if it has been dropped or has fallen into water or has been damaged.
- The scales have a measuring range up to 180 kg, 396 lb/28.8 st. Do not overload the scales.
- Place the scales on a firm, level base. Soft, uneven surfaces are unsuitable for the measuring process and lead to incorrect results.
- Do not use the scale in a location where extremely high temperatures or humidity may occur.

1.2 Important

- Keep the scale away from water.
- Handle the scale with care. The scale must not be subjected to impacts or vibration. Do not drop the scale.
- Do not attempt to dismantle the scale, otherwise the warranty will be invalidated. This unit contains no user-maintainable or user-replaceable parts apart from the battery.
- Do not attempt to repair the device yourself in the event of a malfunction since this will invalidate the warranty. Repairs should only be carried out by authorised service centres.
- If the scale is not going to be used for a long period, please remove the batteries.

2.1 Scope of supply and packaging

Please check first of all that the unit is complete.
The following parts are included as standard:

- 1 **MEDISANA PST** Personal scales
- 1 Type CR2032 3V lithium coin cell battery
- 4 non-slip stickers
- 1 Instructions for use

The packaging can be reused or recycled. Please dispose properly of any packaging material no longer required. If you notice any transport damage during unpacking, please contact your dealer without delay.

WARNING

To avoid risk of suffocation, keep packaging film away from children!

2.2 Performance features

- Top quality safety glass
- Top quality sensor technology
- Measuring range 1 kg to 180 kg or 2.2 lb to 396 lb/28.8 st
- 100 g, 0.2 lb or 1/4 lb graduation
- Switchable between kg, lb and st
- Direct measurement of body weight
- Easy to use
- Clear LCD display
- Includes 1 type CR2032 3V battery

2.3 Prerequisites for correct measuring results

Reliable measurements can only be achieved under certain prerequisites:

- Always weigh yourself under the same conditions.
- Always weigh yourself at the same time of day.
- Stand up straight and stand still.
- Weigh yourself preferably in the morning, after showering or bathing, and with dry feet.

It is possible to obtain a reliable indication of body weight development from the results if measurements are carried out consistently.

3.1 Insert/replace batteries

Before you start up your digital scale, insert the provided cell battery (type CR2032) into the device. To do this, open the battery compartment **1** beneath the device and insert the battery. Please note polarity (positive pole upwards; facing towards the lid). Close the battery compartment again. Replace the battery if the battery change symbol "Lo" appears in the display **3** or if nothing appears in the display when the device is switched on.

BATTERY SAFETY INFORMATION

- **Keep away from children.**
- **Do not short-circuit.**
- **Do not recharge.**
- **Do not throw into fire.**
- **Do not dispose of old batteries with your household waste, but dispose of them at a battery collection station at a recycling site or in a shop.**

3.2 Non-Slip Stickers

4 non-slip stickers are included with the **MEDISANA** personal scales **PST**. For your safety, these should be positioned on the glass plate. Remove the protective backing and stick them on the clean, completely dry glass plate.

3.3 Weighing

Place the scales on a hard, level base. This ensures accurate measuring results and is the requirement for safe weighing without risk of injury.

1. Switch the scales on by briefly tapping the standing area **4** with a foot.
2. "8888" will appear briefly on the display **3** followed by "0.0 kg". The personal scales are ready to measure.
3. Stand on the scales and remain still. "kg" flashes on the display whilst your weight is being determined. "kg" will appear continuously when the measuring result has stabilised. Your measured body weight, e.g. **65.1** (65 kg 100 g) will be displayed for approx. 10 seconds. The scales then switch off automatically. Step down from the scales.

**NOTE**

The scales will switch off automatically to save batteries if you do not step onto the scales within 10 seconds.

3.4 Change the Unit of Weight

The measuring unit of the **MEDISANA** personal scales **PST** may be changed over. Turn the scales over and select the desired measuring unit KG, LB or ST using the black changeover switch **2**.

3.5 Error messages

| | |
|-------------|---|
| Lo | The unit switches off after approx. 2 seconds. Replace the battery. |
| FFFF | The scale is overloaded. |
| ERRO | The weight sensor is faulty. The measuring results displayed are too low. |
| ERR1 | The weight sensor is faulty. The measuring results displayed are too high. |

4.1 Cleaning and maintenance

- Remove the battery before cleaning the unit.
- Never use aggressive cleaning agents or stiff brushes.
- Clean the scale using a soft and lightly moistened cloth. Never use abrasive cleaning agents or alcohol.
- Do not allow any water to get into the unit. Do not use the unit again until it is completely dry.
- Remove the battery from the unit if it will be unused for a longish period. Otherwise there is a risk that the batteries might leak.
- Do not expose the unit to direct sunlight; protect it against dirt and moisture.

4.2 Disposal

This product must not be disposed together with the domestic waste.

All users are obliged to hand in all electrical or electronic devices, regardless of whether or not they contain toxic substances, at a municipal or commercial collection point so that they can be disposed of in an environmentally acceptable manner.

Please remove the battery before disposing of the equipment. Do not dispose of old batteries with your household waste, but dispose of them at a battery collection station at a recycling site or in a shop.

Consult your municipal authority or your dealer for information about disposal.

4.3 Technical data

| | |
|------------------------|--|
| Name and model | : MEDISANA PST Personal scales |
| Dimensions (L x W x H) | : approx. 30 cm x 30 cm x 1.85 cm |
| Weight | : approx. 1.97 kg |
| Voltage supply | : 3 V = , 1 x 3 V lithium coin cell battery type CR2032 |
| Measuring range | : 1 kg to 180 kg or 2.2 lbs to 396 lbs/28.8 st |
| Graduation | : 100 g, 0.2 lb or 1/4 lb |
| Autom. switch-off | : after approx. 10 sec. |
| Item number | : 40420 |
| EAN number | : 4015588 40420 7 |

**We reserve the right to make technical and design changes
in the course of continuous product improvement.**

5.1 Warranty and repair terms

Please contact your dealer or the service centre in case of a claim under the warranty. If you have to send in the unit, please enclose a copy of your receipt and state what the defect is.

The following warranty terms apply:

1. The warranty period for **MEDISANA** products is three years from date of purchase. In case of a warranty claim, the date of purchase has to be proven by means of the sales receipt or invoice.
2. Defects in material or workmanship will be removed free of charge within the warranty period.
3. Repairs under warranty do not extend the warranty period either for the unit or for the replacement parts.
4. The following is excluded under the warranty:
 - a. All damage which has arisen due to improper treatment, e.g. non-observance of the user instructions.
 - b. All damage which is due to repairs or tampering by the customer or unauthorised third parties.
 - c. Damage which has arisen during transport from the manufacturer to the consumer or during transport to the service centre.
 - d. Accessories which are subject to normal wear and tear.
5. Liability for direct or indirect consequential losses caused by the unit are excluded even if the damage to the unit is accepted as a warranty claim.

MEDISANA AG
Itterpark 7-9
D-40724 Hilden
Tel.: +49 (0) 2103 / 2007-60
Fax: +49 (0) 2103 / 2007-626
eMail: info@medisana.de
Internet: www.medisana.de

The service centre address is shown on the last page.

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D&L Marketing Oy
ALERE Electronics AB