MEDISANA®

Shiatsu-Massagegerät SKM

Shiatsu Massager SKM

Appareil de massage Shiatsu SKM

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Aparato de masaje Shiatsu SKM

Aparelho de massagem Shiatsu SKM

Shiatsu-massageapparaat SKM

Shiatsu-hierontalaite SKM

Shiatsu-massageapparat SKM

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1.1 Thank you very much

Thank you for your confidence in us and congratulations on your purchase!
With the Shiatsu-massage unit **SKM** you acquired a quality product from **MEDISANA**.

To ensure the best results and long-term satisfaction with your **MEDISANA SKM** massager, we recommend that you read the following operating and maintenance instructions carefully.

Read the instruction manual carefully before using this device, especially the safety instructions, and keep the instruction manual for future use.

Should you give this device to another person, it is vital that you also pass on these instructions for use.

1.2 Important

- Before connecting the device to your power supply, please ensure that the supply voltage stated on the rating plate is compatible with your mains supply.
- Only use this unit according to the intended purpose and as directed in the instructions for use. The device is not intended for commercial use.
- Never use the unit outdoors.
- Do not sit or stand on the unit; it is not designed to support your full body weight.
- Never cover up the ventilation holes on the bottom of the device.
- Do not put any foreign objects into any of the openings. Avoid contact with pointed or sharp objects. The warranty is null and void if it is not used for the intended purpose.
- If the unit is to be used by or on people who are handicapped in any way or in the presence of children, this must always be done under supervision. The unit is not designed to be used as a toy!

WARNING

A single treatment session should not last for more than 10 - 15 minutes.

- Do not attempt to repair the unit yourself in the event of any defects. Not only does this render the warranty null and void but you can also run the risk of serious hazards (fire, electric shock, injury). Only have repairs carried out by authorised service points.
- If the power supply cord is damaged, it must be replaced by the manufacturer, its service agent or similarly qualified persons in order to avoid a hazard.

Note the following in particular concerning electrical equipment:

- Do not reach for the unit if it has fallen into water. First pull out the mains plug at once.
- Do not use, store or keep the unit in the bathroom. Do not store or keep any electrical equipment in places where it could fall into the bathtub or washbasin.
- Do not use the unit under the bedclothes or under cushions.

1.2

Important

- Keep the power supply cord away from hot surfaces.
- Never use the unit in rooms in which aerosols (sprays) are being used or where pure oxygen is being administered.
- Never use the unit if the power supply cord has been damaged, if the unit does not work perfectly, if it has fallen into water or has been damaged.
- Never carry, pull or rotate the unit by the power supply cord and never pinch the supply cord between other objects.
- Only plug in or unplug the optional remote control when the unit is not connected to the mains outlet.
- Do not allow the unit to become overheated due to excessively long use.
- Unplug the unit from the mains outlet when it is switched off and no longer in use.
- Only pull the plug out of the socket when you have dry hands.

1.3

Important recommendations concerning your health

Please remember that home massage treatment cannot replace medical treatment of the kind that is necessary in the case of serious complaints.

Do not use the unit directly on bones or joints, the head, the front of the neck, or in the genital or kidney areas.

You should refrain from using the **SKM** Shiatsu massager, or should consult your doctor before using it, if

- you are pregnant,
- you have a cardiac pacemaker, artificial joints or electronic implants,
- you suffer from one or more of the following illnesses and complaints:
Circulation disorders, varicose veins, open wounds, bruises, broken skin, inflammation of the veins.
- If you experience pains in the muscles and joints that persist after a long period of time, please inform your doctor. Persistent pains can be symptomatic of serious illness.
- Should you experience pain or discomfort during treatment, stop treatment immediately.
- If you experience unexplained pain or swelling in the legs or feet, or suffer from a muscle injury, you should consult your doctor before using the unit. Stop the massage if you experience pain during or swelling after treatment.

2.1

Scope of delivery

Please check first that the unit is complete.

The package contains:

- 1 **SKM** Shiatsu Massager
- 1 Operating instructions

Accessory (optionally available, not included in set):
remote control for **SKM** Shiatsu Art.-No. 88856.

2.2 Useful information about Shiatsu massage

Shiatsu is a traditional pressure massage technique from the Far East that can be used to stimulate and relax you, from your neck right down to the soles of your feet.

Shiatsu massage predominantly uses finger pressure to relieve tenseness and have a harmonising effect on body and mind. Pressure can be exerted statically, dynamically, on individual points or on wider areas.

The **MEDISANA SKM** Shiatsu massager features two rotating massage heads equipped with an additional vibration function. The combination of vibration and pressure helps loosen up tense muscles. The pressure and vibration massage functions can be activated individually or in combination.

3.1 Getting started

Before plugging in the power plug, pay attention that the function levers for rotation **4** and vibration **5** are fixed in the **0** – position (middle position). Then plug the power plug **2** into the designated socket (230 V).

3.2 Switching on and operating the unit

On the left-hand side of the unit you will find two switches **4** **5** that control the rotation and vibration functions. You can slide both switches to the left and the right of the central position. When in the central position, the respective function is switched off. Both functions can be activated separately or in combination.

Setting rotation

If you slide the rotation selector switch **4** to the left, the two massage heads will rotate forwards. If the switch is moved to the right position, the heads rotate backwards.

Setting vibration

Slide the selector switch **5** to the left to cause the massage heads to vibrate gently. Moving the switch to the right of the central position will give you a stronger vibration.

3.3 Massage tips

You can perform a massage either on yourself or on a partner. Switch on the unit and place it on the part of the body to be massaged. Do not exert any pressure, just gently apply the massage heads to the part of the body being treated. The rhythmic massaging motion and the vibration of the massage heads will help loosen up the muscles. Should you find the massage pressure too strong, you can cover the massage heads with a towel or turn down the massage intensity. To massage the neck region, arms or lower leg, place the unit on a firm and stable surface, cover the massage heads with a towel and place the part of the body to be treated directly on or between the massage heads.

- Treat the parts of the body you wish to massage for 1-3 minutes each.
- A single treatment session should not last for more than 10 - 15 minutes.
- Massage treatment should always be a pleasant experience. If you experience any pain or discomfort during massage, please stop treatment and consult your doctor.

3.4 Switching off the unit

Switch the unit off by returning the selector switches **4** **5** to their original central position. When you have finished using the unit, please unplug it from the mains outlet. After completing a massage session, wait at least 15 minutes before switching on the unit again. Should the unit overheat during use, it will switch off automatically.

3.5 Using the remote control (optional)

The **SKM** Shiatsu massager can be controlled either using the selector switches on the unit or with a remote control unit optionally available from **MEDISANA**. To use the remote control, plug it into the designated socket **3** on the unit, making sure that the plug's catch spring is pointing upwards. Now plug the **SKM** Shiatsu massager into a suitable mains outlet (230 V). A lightning symbol will start to flash on the remote control's display. The remote control has 5 buttons for controlling the various functions, with corresponding symbols being shown on the display.

Button 1 Standby

The flashing lightning symbol on the display indicates that the remote control is connected to the power. After a short press of button 1, the symbol will stop flashing and remain constantly visible. The massager can now be operated using the remote control.

Button 2 Rotation

Pressing this button enables you to set the rotation of the massage heads to forwards or, with a further press of the button, to backwards. Pressing the button one more time switches off the rotation function.

Button 3 Vibration

This button sets vibration intensity. You have the choice of three settings – low – medium – high. Each press of the button selects the next-highest vibration setting. A fourth press of the button switches off massage head vibration.

Button 4 Duration

This button enables you to set the duration of the massage to 5, 10 or 15 minutes. After each press of the button, the set massage duration will appear on the display. You can make your selection by repeatedly pressing the button. Once the duration has been set, the timer will count down to zero, after which the unit will switch off automatically. If you do not wish to set a fixed time or want to stop the massage before the timer has reached zero, press the button one more time.

Button 5 Programme selection

With this button you can select one of three preset programmes. Each press of the button selects the next respective programme.

Programme I The massage direction changes every 15 seconds – the vibration intensity changes every 10 seconds.

Programme II The massage direction changes every 15 seconds – the vibration function is deactivated.

Programme III The massage function is deactivated – the vibration intensity changes continuously.

When you have finished using the unit, press the catch spring on the remote control's plug gently downwards and unplug the remote control from the socket. On the underside of the unit you will find a compartment where you can store the remote control after use, with the display facing downwards.

4.1 Cleaning and maintenance

- Before cleaning the unit, make sure it is switched off and disconnected from the mains. Clean the unit once it has cooled down using a soft, damp cloth.
- Do not use any strong detergents for cleaning, e.g. scouring cream or other cleansing agents, as they may damage the surface of the unit.
- When cleaning, never immerse the unit in water and make sure that no water gets into the unit.
- Only use the unit again when it is completely dry.
- Store the unit in a cool, dry place.
- Take care when rolling up the mains lead to avoid damaging it.

4.2 Disposal

This product must not be disposed together with the domestic waste.

All users are obliged to hand in all electrical or electronic devices, regardless of whether or not they contain toxic substances, at a municipal or commercial collection point so that they can be disposed of in an environmentally acceptable manner.

Consult your municipal authority or your dealer for information about disposal.

4.3 Specifications

Name and model	:	MEDISANA SKM Shiatsu Massager
Power supply	:	230 Volt~ / 50 Hz
Power consumption	:	40 Watt
Dimensions (L x W x H)	:	approx. 34,5 x 25 x 18 cm
Weight	:	approx. 3,4 kg
Article no.	:	88857
EAN Code	:	40 15588 88857 1

Accessory Remote control

Article no.	:	88856
EAN Code	:	40 15588 88856 4

In accordance with our policy of continual product improvement, we reserve the right to technical and design changes without notice.

5.1 Warranty and repair terms

Please contact your dealer or the service centre in case of a claim under the warranty. If you have to send in the unit, please enclose a copy of your receipt and state what the defect is.

The following warranty terms apply:

1. The warranty period for **MEDISANA** products is three years from date of purchase. In case of a warranty claim, the date of purchase has to be proven by means of the sales receipt or invoice.
2. Defects in material or workmanship will be removed free of charge within the warranty period.
3. Repairs under warranty do not extend the warranty period either for the unit or for the replacement parts.
4. The following is excluded under the warranty:
 - a. All damage which has arisen due to improper treatment, e.g. non-observance of the user instructions.
 - b. All damage which is due to repairs or tampering by the customer or unauthorised third parties.
 - c. Damage which has arisen during transport from the manufacturer to the consumer or during transport to the service centre.
 - d. Accessories which are subject to normal wear and tear.
5. Liability for direct or indirect consequential losses caused by the unit are excluded even if the damage to the unit is accepted as a warranty claim.

MEDISANA AG

Itterpark 7-9

40724 Hilden

Germany

Tel.: +49 (0) 2103 / 2007-60

Fax: +49 (0) 2103 / 2007-626

eMail: info@medisana.de

Internet: www.medisana.de

The service centre address is shown on the last page.

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D&L Marketing Oy
ALERE Electronics AB