MEDISANA®

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GB

- 1 Rotary switch with 4 positions
 - 0 Off
 - 1 Vibration massage and infrared light radiation; the water temperature is held constant
 - 2 Spray bath with water jet massage the water temperature is heated
 - 3 Spray bath, vibration and water jet massage with infrared light radiation; the water is heated
- 2 Vibrating massage nipples
- 3 Maximum filling height
- 4 S - formed air spray bar
- 5 Dual water jet massage tap
- 6 Depression for acupressure massage inserts
- 7 Infrared light source
- 8 Aroma centre for aromatic substances and bath salts
- 9 Water outlet notches
- 10 Acupressure massage inserts

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Thank you very much

Thank you for your confidence in us and congratulations!

With the **MEDISANA** foot spray bath **WBM** you have purchased a high-quality product which will contribute to your wellness for a long time if correctly used and serviced.

In order to achieve the desired effect with your foot spray bath in the long term, we recommend that you read the following information on its use and maintenance carefully.

1.1 Important

Please Please read the following safety information carefully before using the device and keep these instructions in a safe place for later use.

- Before connecting the device to your power supply, please ensure that the supply voltage stated on the rating plate is compatible with your mains supply.
- Do not immerse the device in water or other fluids.
- If, nonetheless, liquid should penetrate the device, unplug the mains cable immediately.
- Do not use the device if the cable or socket is damaged and keep the cable away from hot surfaces.
- Only use the device with water and appropriate bath oils and not with other liquids.
- Turn the device off after use and pull the plug out of the socket.
- Do not stand up inside the device as it is not designed to take the burden of full body weight.
- Place the device on a firm floor area when in use so that it cannot tip over.
- Touch the switch and plug only with dry hands.
- Do not use the cable to pull or carry the device.
- Pay special attention to correct use if the device is used by or in the vicinity of children or handicapped persons.
- The device is not determined for commercial use.
- Do not leave the device unattended when it is switched on.
- Only use the device for its intended use as in the instructions for use.
- Do not use the device outdoors or in damp rooms.
- Never cover up the ventilation holes on the bottom of the device
- If the device malfunctions please do not repair it yourself. Do not use the device further. Ask your specialist dealer and only have repairs carried-out by authorised service companies.
- Should you experience medical complaints during use talk to your doctor. Should this be the case, stop using the device immediately.
- If you have medical problems with your feet, vein diseases or diabetes, talk to your doctor before using the device.

1.2 Recommendations for use

Duration

- Generally speaking you can take about 10 to 15 minutes for a footbath with a massage, once or twice daily.
- However you should ensure that the device has completely cooled off before it is used again!

General health recommendations

- Everybody responds differently to a massage. It varies from person to person whether a gentle or an intensive massage is found to be pleasant. After using the **MEDISANA WBM** foot spray bath a few times you will discover your own personal optimum application form to achieve the desired relaxation.
- Should you have any health questions regarding the application, your doctor will certainly be able to help you further.
- Should you have inexplicable pain or swellings in your legs or feet or have injured a muscle, please consult your doctor before use.
- Stop a massage if you experience pain during the application or swellings crop up after use.
- We advise consultation with your doctor if you are pregnant.
- If you suffer from Diabetes or diseases of the veins, please consult your doctor before use.

2.1 Pedicure using the WBM foot spray bath

Well cared for feet not only look good, they are generally speaking also healthy. Just a few minutes' attention per day suffice for highly stressed feet to feel good and cope with the daily chores without pain. One important part of pedicure is a relaxing footbath.

The **MEDISANA WBM** foot spray bath offers the best pre-requisites for this. It makes a powerful foot massage with or without water possible, an intensive reflex zone massage and infrared light heat radiation too. If required, the device can be used for the massage or spray function only or it can also be combined with the infrared radiation. Whatever you select, the varied applications help you to achieve an effect for body and soul which most people find pleasant and relaxing.

Massage

Equipped with vibrating nipples, an air spray bar, massage jets and exchangeable acupressure inserts, the **MEDISANA WBM** foot spray bath enables a foot massage from the sole to the foot instep.

The foot massage supplies the tissue layers better with blood and can thus unblock veins and lymphoid vessels. Tense muscles can be loosened and painful joints soothed.

The metabolism in the cells is stimulated.

Infraredradiation

The infrared radiation turned on in positions **1** and **3** penetrates the soles of the feet pleasantly and can make a far-reaching heat therapy possible.

Heattreatment

Das **MEDISANA WBM** foot spray bath heats the water in the device to a pleasant temperature and ensures that the temperature remains constant during use.

Aromatherapy

A depression is embedded in the **MEDISANA WBM** footbath for the addition of aromatic substances and bath salts.

The added substances mix slowly with the water and intensify the relaxing effect of the footbath.

3.1 Scope of delivery

Please check first whether the device is complete.
The following belong to the scope of delivery:

- The **MEDISANA** foot spray bath **WBM**
- 4 Massage inserts
- Instructions for use

3.2 Footbath and massage with water

Before filling the device with water, make sure that no mains leads are connected, i.e. the electricity lead must not be plugged-in!
When the foot spray bath is securely placed on the floor (ensuring that it cannot tip over), and has been filled with water, insert the mains plug into the power supply socket.

3.3 Foot massage without water

You can also use the device for a foot massage without water - for this mode switch the Rotary switch **1** to position **1** (vibration/massage with infrared radiation).

3.4 Filling with water

The **MEDISANA WBM** foot spray bath can be filled with up to 3.3 litres of water. The maximum filling height is marked **3** in the foot spray bath. You can use cold or warm water with or without bath additives, just as desired.

Never use any foaming bath additives or bath salts.

3.5 Function setting

Sit down comfortably in front of the foot spray bath and select one of the three different functions using the Rotary switch **1**.

Heat function selection

The **MEDISANA WBM** foot spray bath includes an innovative heating function. This technical improvement makes it possible to keep the water at a constant temperature during application or to warm it.

Position 1 : Maintain the water temperature

If you have switched the Rotary switch to **1**, the temperature of the filled water will remain – or it will be heated a little.

Position 2 and 3 : Heating water

If set to position **2** or **3**, a warming function is turned on automatically. Tap water reaches a pleasant temperature for a relaxed footbath after ca. 10 minutes and this is kept constant during application.

If the water in position **2** or **3** is heated too hot for you, turn down to position **1** and wait until the water has cooled down to the temperature you desire. If you add cold water to lower the temperature, please make sure not to exceed the maximum filling height **3**.

ATTENTION

If you are very sensitive to heat, have cardiovascular problems, swellings, inflammations, skin irritation or skin wounds, please consult your doctor before use.

**3.6
Aromatic
substances**

Remove the cap and place the aromatic substances and bath oils in the depression **8** located below. During operation of the device, the additive will slowly be mixed into the water.

**4.1
Cleaning and
maintenance**

- Before cleaning the device, make sure that it is turned off and the mains plug is removed from the socket. Pour the water out over the water outlet notches **9** and clean the device with a cloth when it has cooled off.
- Please do not use any aggressive substances, e.g. scouring chemicals or other washing-up liquids when cleaning. These could damage the surface.
- Store the device in a dry, cool place.
- Wind-up the electricity lead carefully to avoid it parting.

**4.2
Disposal**

Dispose of your old electrical unit in an environmentally acceptable manner when it has reached the end of its service life!

Packaging is reusable or may be recycled; metal parts must be disposed of as scrap metal; plastics, electrical and electronic components must be disposed of as electroscrap. Consult your municipal authority or your dealer for information about disposal.

**4.3
Specifications**

Name und model	: MEDISANA foot spray bath WBM
Dimensions L x W x H	: 47 x 40 x 23 cm
Weight	: 3.2 kg
Current supply	: 230 Volt / 50 Hz
Power consumption	: 490 Watts
Article no.	: 88362
EAN-Number	: 4015588883620

In accordance with our policy of continual product improvement, we reserve the right to technical and design changes without notice.

5.1 Warranty and repair terms

Please contact your dealer or the service centre in case of a claim under the warranty. If you have to send in the unit, please enclose a copy of your receipt and state what the defect is.

The following warranty terms apply:

- 1.** The warranty period for **MEDISANA** products is three years from date of purchase. In case of a warranty claim, the date of purchase has to be proven by means of the sales receipt or invoice.
- 2.** Defects in material or workmanship will be removed free of charge within the warranty period.
- 3.** Repairs under warranty do not extend the warranty period either for the unit or for the replacement parts.
- 4.** The following is excluded under the warranty:
 - a.** All damage which has arisen due to improper treatment, e.g. non-observance of the user instructions.
 - b.** All damage which is due to repairs or tampering by the customer or unauthorised third parties.
 - c.** Damage which has arisen during transport from the manufacturer to the consumer or during transport to the service centre.
 - d.** Accessories which are subject to normal wear and tear.
- 5.** Liability for direct or indirect consequential losses caused by the unit are excluded even if the damage to the unit is accepted as a warranty claim.

Please do not hesitate to contact us if you have questions about this or any other product.

MEDISANA AG

Itterpark 7-9

D-40724 Hilden

Tel.: +49 (0) 2103 / 2007-60

Fax: +49 (0) 2103 / 2007-626

eMail: info@medisana.de

Internet: www.medisana.de

The service centre address is shown on the last page.

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D&L Marketing Oy
ALERE Electronics AB